



Calibration Manual

Oscilla® SM960-C

Table of content:

1. GAIN CONSTANT adjustment (mode1).....	4
2. INPUT adjustment (mode 2).....	5
3. OUTPUT adjustment (mode 3).....	6
4. LEVEL OFFSETS adjustment (mode 4)	8
5. SENSOR CORRECTIONS (mode 5).....	11

Due to tolerances of the transducers and the hardware, it is necessary to adjust the audiometer's output signal to the correct level. The adjustment is done with the telephones, bone conductor, insert phones and freefield speaker the audiometer is used with. If one of the transducers is replaced, the audiometer should be adjusted again, with the new transducer. Adjustment of tone should be done through a bandpass filter, which is set to the frequency that is adjusted. Adjustment of masking is done without bandpass filter. Before entering the calibration mode, connect a tone generator to the TAPE input, select the TAPE input on channel 1 and turn up the tone generator's output level to the point just where the 0 dB bar in the VU display turns on. The tone generator's frequency should be 1000 Hz, sine wave.

To enter the calibration mode, do the following:

Turn off the audiometer.

Press and hold down the keys 'ATTENUATOR DOWN', 'INTR.' and 'FUNC.' on channel 1.


Turn on the power and wait for approx. 2 second. Release the three keys again.

The 'response' lamps will flash 5 times, and stay on. The VU display will show 'CALIBRATION', and then 'GAIN CONSTANT'. The audiometer is now in calibration mode.

There are five calibration modes:



- 1: GAIN CONSTANT
- 2: INPUT
- 3: OUTPUT
- 4: LEVEL OFFSETS (RETSPL)
- 5: SENSOR CORRECTIONS (Default, not in use)

Use the channel 2 attenuator to choose between the four modes. The calibration mode number is shown in the right display. Calibration must be done in sequence 1, 2, 3. Mode 4 is not really a calibration, but sets and displays the current hearing level standard for the different transducers.




1. GAIN CONSTANT adjustment (mode1)

This adjustment is to ensure that the attenuators attenuation is accurate at both large and small attenuations. The adjustment is done at 1000 Hz through bandpass filter. Use the frequency keys to choose between 20, 40, 60, 80 and 100 dB. The adjustment is done with the channel 1 attenuator. A number in the range 0 to 199 is shown in the VU display to show where the adjustment is standing. There is only one adjustment for each channel, which will increase or decrease the difference between the 5 sound levels.



Adjustment

Read the sound level at the five different levels. Adjust so that the difference is 20 dB. The calibration number indicates the difference between the 20 dB steps, so if the difference is too small, adjust up, and if too large adjust down. Repeat the reading of the 5 sound levels and the adjusting, to get as close as possible to the 20 dB difference. Select the opposite channel with the channel 1 output selector, and repeat the adjustment with the other telephone.



2. INPUT adjustment (mode 2)

This is used to adjust the five different input signal levels in relation to each other. Choose the signal type with the channel 1 input selector, and do the adjustment with the channel 1 attenuator:



TONE:

Should not be adjusted, normally standing at 15 dB. If needed the TONE may be increased, but notice that this will affect all places where the tone is used, and therefore requires that the other signal levels (masking and tape) and OUTPUT mode, is readjusted.



Adjust so the signal from TAPE/CD has the same level as the tone. Ensure that the sensitivity is set to 0 dB deflection on the VU-meter first, as mentioned in the first section, by using the Vol. buttons on the front. The input signal in the tape plug, could be a CD-player with the calibration track (1000 Hz sinus tone) playing.

Note: Some transducer manufacturers state a specific sound level offset in reference to the tone, for speech audiometry. For example, for the Eartone 3A and 5A insert phones it is stated that the level should be 12.5 dB above the tone, in a 2cc coupler. This may be obtained by adjusting the TAPE/CD input 12.5 dB up.



NB (Narrow Band) adjustment:

Adjust to 3 dB higher level than the tone at 1000 Hz. Adjust without bandpass filter.






SN (Speech Noise) adjustment:

Adjust to 3 dB higher level than the tone at 1000 Hz. Adjust without bandpass filter.



WN (White Noise) adjustment:

Adjust to 3 dB higher level than the tone at 1000 Hz. Adjust without bandpass filter.



3. OUTPUT adjustment (mode 3)



Used to calibrate the levels of telephones, bone conductor, insert phones and free field speaker. Choose the output to calibrate with the output selectors. The right, left, bone and mono insert outputs are selected on channel 2, and free field is selected on channel 1.

The right and left outputs has two different sets of calibration. One for normal phones, usually TDH-39, and one for inserts, like the Eartone 3A or 5A. To select between these two sets, press the TALK THROUGH button.

In addition to that,

Theory of operation

The output calibration is made so that all frequencies should be calibrated to the same sound pressure level. The audiometer will subtract a hearing level value during calibration, so when it returns to normal operation, the output sound levels will be the correct hearing levels. These hearing levels are stored in memory, and can be displayed and altered in calibration mode 4. Normally they should not be changed, but if the ISO standard is changed, or another set of hearing levels is desired for some reason (different countries may have different standards), they can be changed. If they are changed, it must be done before the output calibration in mode 3 is performed, otherwise the change will have no effect.



Adjustment of the telephones


Tone

Select 125 Hz and left channel. Set the bandpass filter to 125 Hz. Place the left telephone on the microphone. Adjust the telephone's position on the microphone to get the biggest deflection. This will ensure that the telephone has a tight connection to the microphone.



Adjust sound level:

Adjustment of sound level is done at 80 dB, as default. However, if another calibration level is desired, it may be changed by holding down the SHIFT button while pressing the channel 1 attenuator keys. The actual calibration is done with the attenuator keys, without holding SHIFT. When calibrating, a dB udstyrnumber is shown in the VU display. This is the calibration offset for the selected frequency and transducer, and is changed in 0.5 dB steps. Because the hearing levels set in mode 4 are subtracted here in calibration mode, the output levels should be set to the same level (e.g. 80 dB) at all frequencies. This makes it faster and more practical to calibrate, because you don't have to change the sensitivity of the calibration equipment for each frequency. Use the frequency keys to change frequency, and the output 2 output selector to change between left and right channel.

The audiometer has RETSPL values for two telephone types, the TDH39 and the DD45. To select which one is used, press the "ATT. LOCK" button. This will toggle between the two types, and display the type selected.



80 dB SPL



4. LEVEL OFFSETS adjustment (mode 4)

RETSPL table for TDH39 telephone (ISO-389-1995)

Hz	dB
125	45
250	25.5
500	11.5
750	7.5
1000	7
1500	6.5
2000	9
3000	10
4000	9.5
6000	15.5
8000	13

RETSPL table for DD45 telephone

Hz	dB
125	45.5
250	28
500	14
750	7.5
1000	7
1500	9
2000	9.5
3000	11
4000	12.5
6000	18.5
8000	17

Narrow band masking

The masking level should be 3 dB higher than the tone's, at all frequencies. Adjust the masking level without bandpass filter. As the masking level is changed automatically when the tone is adjusted, it is normally not necessary to adjust, but due to the telephones' frequency response at high frequencies, adjustment should be done at least at 4000, 6000 and 8000 Hz. Choose between tone and noise with the channel 2 input selector.

Adjustment

Read the tone's sound level. Select noise with the channel 2 input selector. Adjust, with the channel 1 attenuator, so that the average level is 3 dB higher than the tone's. Repeat the adjustment of tone and masking with the right telephone. Use the channel 2 output selector to choose right channel.

Adjustment of the stereo insert phones

A set of stereo insert phones may be connected to the same outputs as the normal phones. To calibrate them, press the TALK THROUGH button. This toggles between the phones and the inserts. A brief text (PHON or INST) is showed in the frequency display, to indicate which are selected. In addition, when the TALK THROUGH lamp is on, the inserts are selected. The calibration procedure is the same as described in "Adjustment of the telephones".

The hearing level offsets for inserts are different from those used with the phones. The table below shows the hearing levels for Eartone 3A or 5A insert phones when using the DB-138 2cc coupler, given in the ISO 389 standard. These levels are stored in mode 4, as default. They are used for both the stereo insert set and for calibration of the mono insert output. If another coupler is used for calibration, the values should be corrected to fit the table provided with the insert phones.

Hearing level table for Eartone insert phones

Hz	dB
125	26.0
250	14.0
500	5.5
750	2.0
1000	0.0
1500	2.0
2000	3.0
3000	3.5
4000	5.5
6000	2.0
8000	0.0

Adjustment of the bone conductor

Adjustment of the bone conductor is done at 60 dB, as default. As for the telephones, this default level can be changed by holding SHIFT while pressing the attenuator keys. The adjustment should be done through a bandpass filter. At 8000 Hz, no adjustment is made, as this frequency is unused with the bone conductor.

Adjustment

Place the bone conductor on the vibration sensor. Choose 'BONE' with the channel 2 output selector. Adjust to 60 dB at each frequency. As for the telephones, the audiometer will subtract the hearing levels stored in mode 4 to give the same sound pressure at all frequencies while calibrating.

Hearing level table for bone conductor (RETFL)

Hz	dB
125	82.5
250	67.5
500	58.0
750	48.5
1000	42.5
1500	36.5
2000	31.0
3000	30.0
4000	35.5
6000	40.0

Narrow band masking

Perform the same adjustment for narrow band masking levels, as for the telephones.

Adjustment of the mono insert telephone

The mono insert phone should be adjusted to the same levels as the normal telephones. The procedure is the same as for the telephones. Use the output 2 output selector to select the mono insert phone. Calibration of insert phones requires a special acoustic coupler, for example the Brüel & Kjær DB-138 2cc coupler.

Adjustment of the free field speaker

The calibration of tone to the free field speaker is the same on the audiometer, as for the telephones. The only difference is that the audiometer automatically enables the warble tone when free field is selected. This is to give some compensation for the reflections from the test room, so that the physical speaker and microphone position doesn't affect the sound level too much. Refer to the free field calibration manual for a more detailed description of this. The warble tone may be enabled or disabled by means of the WARBLE button, while in calibration mode 3.

Hearing level table for free field

Hz	dB
125	22
250	11
500	4
750	3
1000	2
1500	0
2000	-1.5
3000	-4
4000	-6.5
6000	2.5
8000	11.5

Hearing level offsets (mode 4)

If, for some reason, it is desired to change the hearing level offsets for the different transducers, it is done in this calibration mode. When in this mode, no sound is output, it is only a mode to change and view the hearing level offsets used to subtract when calibration in mode 3 is performed. The changes are done using the same procedure as when calibrating. Use the frequency keys to change frequency, the output 1 and 2 selectors to choose transducer, and attenuator 1 to adjust the selected hearing level offset. It is recommended to check the numbers in mode 4 to make sure they correspond to the hearing level standard in use, before performing a mode 3 calibration. The values shown when selecting insert with the channel 2 output selector, are used for stereo insert phones too.

5. SENSOR CORRECTIONS (mode 5)

It is possible to add corrections for the calibration microphone(Artificial Ear) and Bone mastoid. But as default the function is disabled.

Ending calibration

All adjustments are stored while they are made, in a non-volatile memory. To exit adjustment mode, press the 'DATA' key, or simply turn off the power.